

**Figure 1**

<b><u>Allosteric Effector</u></b>	<b><u>Structure or Name</u></b>
<b>PCHA-DPG</b>	Penta-cyclohexylammonium 2,3-diphosphoglyceric acid
<b>5Na-DPG</b>	Penta-sodium 2,3-diphosphoglyceric acid
<b>IHP</b>	Inositol hexaphosphate
<b>CHA</b>	Cyclohexylammonium
<b>CHA-IHP</b>	Cyclohexylamine added to IHP to give a solution with a pH = 7.1-7.4

**Figure 2**

Effector	P <sub>50</sub> CONTROL WB mmHg	P <sub>50</sub> EFF: WB mmHg	CONC. EFF mM	CONC EFF:WB mM	OSMOL. EFF mOsM	pH EFF.	pH EFF:WB or EFF: <i>fHb</i>	Volume Ratio EFF:WB
HBS+					310		7.22	
<i>fHb</i> in HBS+							7.22	
<b>PCHA-DPG</b>								
WB	37	55.5	30	22	205	7.8	7.42	1:0.375
WB	25	36	30	22	221		7.33	1:0.375
WB	37	37	30	22	313		7.12	1:0.375
WB	37	37	30	12	317		7.02	1:1.5
WB in Bis-Tris	37	38.2	30	22	341		7.1	1:0.375
WB	37	36	30	12	310		6.8	1:1.5
<b>5Na-DPG</b>								
WB	37	38.2	30	22	163		7.4	1:0.375
WB	37	39.5	30	22	321		7.1	1:0.375
<b>IHP</b>								
WB	37	38.2	30	22	185		7.3	1:0.375
<i>fHb</i>	16	40.5					7.19	0.25 $\mu$ M EFF
<b>CHA</b>								
WB	26.8	28.5	30	22	220		6.23	1:0.375
WB	26.8	26.8	30	22	245		6.75	1:0.375
<b>CHA-IHP</b>								
WB	26.8	42	30	22	220		6.36	1:0.375
WB	24.7	58.2	25	14.3	171		6.93	1:0.35
<i>fHb</i>	16	53					7.17	0.25 $\mu$ M EFF

*fHb* = free hemoglobin; WB = whole blood; EFF = effector.

**Figure 3**

<u>Sample</u>	<u>Observed O<sub>2</sub> P<sub>50</sub> (torr)</u>
Human whole blood	9.3 (pH 7.47)
Washed Goldfish blood cells	20.0 (pH 7.52)
Human fHb	4.7 (pH 7.1)
Goldfish fHb	8.5 (pH 7.1)
Goldfish fHb + 0.25 $\mu$ mol IHP	15.0 (pH 7.1)
Goldfish fHb + 0.5 $\mu$ mol PCHA-DPG	10.3 (pH 7.1)
Goldfish fHb + 0.5 $\mu$ mol ATP	21.0 (pH 7.08)
Goldfish fHb + 0.5 $\mu$ mol GTP	23.0 (pH 7.11)

The data presented in this figure was acquired at 25 C in Bis-Tris buffer at pH 7.2-7.4.


Figure 4

A:	Human whole blood	pH 7.47	$P_{50} = 9.3$
B:	washed goldfish blood	pH 7.52	$P_{50} = 20.0$

All experiments were in Bis-Tris buffer and 25 °C.

Collection methods for blood from up to 20 goldfish requires higher amounts of anticoagulant than from a single human source. From published protocols, a special washing at 4 °C and additional steps to remove nucleic acids from lysed red cells are required.

Obviously the O<sub>2</sub> dissociation for Fish hemoglobin entrapped within a red cell is optimal at a lower temperature than Human .






Figure 5

The experiment in figure 4 were recorded at a pH optimal for humans. The pH for fish studies is generally lower, pH 7.1 versus 7.4. Previously we have reproduced data showing the strong pH dependence for human free hemoglobin at different pH. These experiments (fig. 10 and 11) are at pH 7.1. Because of variations with isolating fish red cells and the temperature coefficient on pH illustrate these effects. Previously I had acquired some familiarity with HEPES Buffered and adjusting for its temperature coefficient, they did not translate to Bis-Tris buffers.

A:	human free hemoglobin	pH 7.1	P <sub>50</sub> = 4.7
B:	goldfish free hemoglobin	pH 7.1	P <sub>50</sub> = 8.5
C:	sample B + 0.25 μmol IHP	pH 7.1	P <sub>50</sub> = 15.0

All experiments were in Bis-Tris buffer and 25 °C.





Figure 6

A:	goldfish free hemoglobin	pH unk*	P <sub>50</sub> = 4.0
B:	sample A + 0.5 μmol PCHA•DPG	pH 7.1	P <sub>50</sub> = 10.3
C:	goldfish free hemoglobin + 0.5 μmol ATP	pH 6.84**	P <sub>50</sub> = 43.0
D:	goldfish free hemoglobin + 0.5 μmol GTP	pH 7.11	P <sub>50</sub> = 23.0
**E:	goldfish free hemoglobin + 0.5 μmol ATP	pH 7.08	P <sub>50</sub> = 21.0

All experiments were in Bis-Tris buffer and 25 °C.

\* pH measured after the addition of PCHA•IHP in cuvette A. All other samples began with a new aliquot of free fish hemoglobin prior to the addition of the allosteric effector.

\*\* After adding the ATP the pH of the sample (Run C) was too low. The pH of the ATP was adjusted to yield an appropriate pH after addition to the sample, Run E



**Fig 7** Oxygen Dissociation Curves of Whole Blood treated with a solution of pentacyclohexylammonium-2,3 diphosphoglyceric acid (PCHA-DPG) and Sodium salt of DPG (PNa-DPG).

**A:** Control (25μL WHOLE BLOOD).

$$P_{50} = 37.0$$


**C:** 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PCHA-DPG. After incubation the system was washed 4X and 15μL RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$$P_{50} = 50.5$$

**E:** 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PNa-DPG.

$$P_{50} = 38.2$$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.



**Fig 6** Oxygen Dissociation Curves of Whole Blood treated with a solution of pentacyclohexylammonium-2,3 diphosphoglyceric acid (PCHA-DPG)

**A:** Control (25 $\mu$ L WHOLE BLOOD).

$$P_{50} = 25.5$$

**F:** 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM PCHA-DPG. After incubation the system was washed 4X and 15 $\mu$ L RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$$P_{50} = 47.3$$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.






Fig9. Oxygen Dissociation Curves of Whole Blood treated with a solution of pentacyclohexylammonium-2,3 diphosphoglyceric acid (PCHA-DPG)

A: Control (25μL WHOLE BLOOD).

$$P_{50} = 25.0$$


B: 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PCHA-DPG. After incubation the system was washed 4X and 15μL RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$$P_{50} = 36.0$$

G: 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PCHA-DPG..

$$P_{50} = 40.0$$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.



**Fig/0** Oxygen Dissociation Curves of Whole Blood treated with a solution of Sodium Salts of DPG and IHP.

- A: Control (25μL WHOLE BLOOD).  
P<sub>50</sub> = 37.0
- C: 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PNa-DPG. After incubation the system was washed 4X and 15μL RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.  
Hypotonic. Osm: 163mOsM  
P<sub>50</sub> = 37.5
- D: 75μL Whole Blood incubated (2-5 min) with 200μL 30mM PNa-DPG. Isotonic. Osm: 321 mOsM  
P<sub>50</sub> = 39.6
- E: 75μL Whole Blood incubated (2-5 min) with 200μL 30mM Na-IHP  
Hypotonic. Osm: 185mOsM  
P<sub>50</sub> = 37.5

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.




Fig 11 Oxygen Dissociation Curves of Whole Blood treated with a solution of Cyclohexylammonium (CHA) and CHA salt of IHP

- A: Control (25 $\mu$ L WHOLE BLOOD).  
 $P_{50} = 26.8$
- B: 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA-IHP.  
 After incubation the system was washed 4X and 15 $\mu$ L RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.  
 $P_{50} = 42.0$
- C: 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA.  
 $P_{50} = 28.5$
- D: 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA.  
 $P_{50} = 26.8$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.

Set 3




Fig 12 Oxygen Dissociation Curves of Whole Blood treated with a solution of Cyclohexylammonium-Inositol Hexaphosphate (CHA-IHP)


D: Control (25 $\mu$ L WHOLE BLOOD).

$P_{50} = 24.7$

C: 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA-IHP. After incubation the system was washed 4X and 15 $\mu$ L RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$P_{50} = 58.2$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.



**Fig/3** Oxygen Dissociation Curves of Whole Blood treated with a solution of Cyclohexylammonium-Inositol Hexaphosphate (CHA-IHP)


**A:** Control (25 $\mu$ L WHOLE BLOOD).

$$P_{50} = 23.5$$

**C:** 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA-IHP. After incubation the system was washed 4X. Whole Blood Cell Pellet was stored for 48 hrs at 4-8°C and 15 $\mu$ L RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$$P_{50} = 50.5$$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.



Fig/4 Oxygen Dissociation Curves of Whole Blood treated with a solution of Cyclohexylammonium-Inositol Hexaphosphate (CHA-IHP)

A: Control (25 $\mu$ L WHOLE BLOOD).

$$P_{50} = 24.8$$

C: 75 $\mu$ L Whole Blood incubated (2-5 min) with 200 $\mu$ L 30mM CHA-IHP. After incubation the system was washed 4X and 15 $\mu$ L RBC were used for measurement of the Hb-O<sub>2</sub> dissociation curve at 37°C.

$$P_{50} = 32.8$$

Incubation Time: 2-5 min at 37°C. All Experiments were conducted with Whole Blood.